

Firmware Update Preparation

In advance of your scheduled firmware update, you can do a few things to better protect freezer contents.

- Locate the rocker ON|OFF switch (#6a below) on the back of the freezer. Flip the switch to the OFF position. Wait for the main display (#3 below) to shut down, then flip the switch back ON. This will not affect any settings or validation performed on the freezer.
- Use a temporary secondary monitoring device to actively read the freezer temperature and alert you if the cabinet set point temperature is not being maintained.

If you have any questions or concerns about your freezers please contact our technical service department at 740.274.7980 for more help.

