

# The Personal Minicentrifuge



## The Personal Minicentrifuge

Take a quick spin with your own Personal Minicentrifuge! The sleek design and compact footprint is the ideal addition to your workbench or portable enough for use in a fume hood. It's your perfect lab partner for quick spin downs of condensate, vortexed or thawed samples. The Personal Minicentrifuge is ideally suited for applications requiring relatively low centrifugal forces, such as microfiltration and cell separation.

Included with your Personal Minicentrifuge unit are a PCR strip rotor, a microcentrifuge tube rotor, a power adapter with 4 exchangeable plugs, and adapters for both 0.5 mL and 0.2 mL tubes.

8030-5000      Personal Minicentrifuge      ~~\$350.00~~

### Specifications:

|               |   |
|---------------|---|
| Max Speed:    | 6,500 rpm $\pm$ 1%  |
| Max RCF:      | 2,500 $\times$ g  |
| Rotors:       | 8 $\times$ 1.5/2.0 mL and 0.5/0.2 mL adapters for microcentrifuge tubes and 4 $\times$ 8-tube strips or 32 individual PCR tubes |
| Noise Level:  | 55 dB(A)  |
| Power Supply: | 100 - 240V, 50Hz - 60Hz   |
| Unit Input:   | 24 VDC 12 W   |
| Weight:       | 2.5 lbs.  |
| Dimensions:   | 5.1" W $\times$ 6.3" D $\times$ 4.4" H  |
| Warranty:     | 2 Years   |

~~\$350.00~~  
\$290



**Quick Start** - The Personal Minicentrifuge is easy to start. Just press down on the clear protective lid to securely lock and start spinning!



**Quick and Easy Rotor Change** - To mount one of the two rotors, simply push it onto the drive shaft until there is an audible click. Remove it by pulling it straight up and off the shaft.



**Quiet Operation** - The quiet and smooth operation does not contribute to laboratory background noise.



**One Touch Stop** - The Personal Minicentrifuge is just as easy to stop. Lightly press the stop button on the right for fast braking and release of the electromagnetic lid lock.