

The Personal Minicentrifuge

Centrifugation at your fingertips

Spin tubes and strips from 0.5 mL to 2 mL

Easy one-handed operation

Take a quick spin with your own Personal Minicentrifuge! The sleek design and compact footprint is the ideal addition to your workbench or portable enough for use in a fume hood. It's your perfect lab partner for quick spin downs of condensate, vortexed or thawed samples. The Personal Minicentrifuge is ideally suited for applications requiring relatively low centrifugal forces, such as microfiltration and cell separation.

The Personal Minicentrifuge unit includes a PCR strip rotor, a microcentrifuge tube rotor, a power adapter with 4 exchangeable plugs, and adapters for both 0.5 mL and 0.2 mL tubes.

8030-5000 Personal Minicentrifuge

~~\$350.00~~

FLASH SALE
ONLY \$279!!



Quick Start - The Personal Minicentrifuge is easy to start. Just press down on the clear protective lid to securely lock and start spinning!



Quick and Easy Rotor Change - To mount one of the two rotors, simply push it onto the drive shaft until there is an audible click. Remove it by pulling it straight up and off the shaft.



One Touch Stop - The Personal Minicentrifuge is just as easy to stop. Lightly press the stop button for fast braking and release of the electromagnetic lid lock.



Quiet Operation - The quiet and smooth operation does not contribute to laboratory background noise.

Specifications:

Max Speed:	6,500 rpm \pm 1%
Max RCF:	2,500 \times g
Rotors:	8 \times 1.5/2.0 mL and 0.5/0.2 mL microcentrifuge tubes adapters and 4 \times 8-tube strips or 32 individual PCR tubes
Noise Level:	55 dB(A)
Power Supply:	100 - 240V, 50Hz - 60Hz
Unit Input:	24 VDC 12 W
Weight:	2.5 lbs.
Dimensions:	5.1" W \times 6.3" D \times 4.4" H
Warranty:	2 Years



8030-5011



8030-5010

A PCR strip rotor, microcentrifuge tube rotor, and tube adapters are included with purchase of a Personal Minicentrifuge. Extras can be purchased separately online.